



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RAC2 Anti-Human Antibody (Clone AT2G11)

单克隆抗体

产品货号: MA01526

第三版

**免疫原:**Recombinant human RAC2(1-189aa) protein purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT2G11

**描述:**Rac2(Ras-related C3 botulinum toxin substrate 2), a small (21 kDa) signaling G protein, belongs to the family of Rho GTPases regulating a variety of biological response pathways that include cell motility, cell division, gene transcription, and cell transformation.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF