

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-ORF73/HHV8 Antibody (Clone AT4C11)

单克隆抗体

产品货号: MA01610

第三版

**免疫原:**Recombinant human ORF73/HHV8 (122-329aa) purified from E. coli

**亚型:**Mouse IgG2a heavy chain and kappa light chain

**克隆:**AT4C11

**描述:**ORF73 protein of human herpesvirus 8(HHV-8) is a Kaposi's sarcoma-associated herpesvirus (KSHV) latent protein which has been suspected of playing a crucial role in modulating viral and cellular gene expression. ORF73 is associated with cellular chromatin and stays on the chromosomes during cell division.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF