

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PKCK2 Anti-Human Antibody (Clone 3G1)

单克隆抗体

产品货号: MA01694

第三版

免疫原: Recombinant human PKCK-2 (1-391aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 3G1

描述: Protein Kinase Casein Kinase 2 (PKCK-2) is a ubiquitous Ser/Thr kinase expressed in all eukaryotes. It phosphorylates a number of proteins involved in various cellular processes. PKCK-2 holoenzyme is catalytically active tetramer, composed of two homologous or identical and constitutively active catalytic () and two identical regulatory () subunits.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF