

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD71 Monoclonal Antibody

单克隆抗体

产品货号: MA01910

第三版

免疫原:Fusion protein of CD71

亚型:IgG2a

克隆:30

描述:CD71,also known as transferrin receptor protein 1 (TfR1),is a transmembrane glycoprotein composed of two disulfide-linked monomers,each of 90 kDa molecular weight. Each monomer binds one holo-transferrin molecule creating an iron-Tf-TfR complex which enters the cell by endocytosis. CD71 is present on actively proliferating cells and is essential for iron transport into proliferating cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF