



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CDK2 Monoclonal Antibody

单克隆抗体

产品货号: MA01913

第三版

免疫原:Fusion protein of CDK2

亚型:IgM

克隆:154

描述:CDK2(Cyclin-dependent kinase 2) is also named as CDKN2 and belongs to the protein kinase superfamily,CMGC Ser/Thr protein kinase family,CDC2/CDKX subfamily.It is involved in the control of the cell cycle; essential for meiosis,but dispensable for mitosis.It has 2 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF