



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CST3 Monoclonal Antibody

单克隆抗体

产品货号: MA02465

第三版

免疫原:Fusion Protein of Cystatin C

亚型:IgG2a

克隆:2E6F3

描述:Cystatin C is a 13-kDa inhibitor of cysteine proteinases which is secreted by all cell types and is completely cleared from the organism through glomerular filtration, shown to be an early and sensitive biomarker of renal dysfunction. A mutation in Cystatin C has been associated with amyloid angiopathy. Expression of this protein in vascular wall smooth muscle cells is severely reduced in both atherosclerotic and aneurysmal aortic lesions, establishing its role in vascular disease. In addition, this protein has been shown to have an antimicrobial function, inhibiting the replication of herpes simplex virus.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF