



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EIF2A Monoclonal Antibody

单克隆抗体

产品货号: MA02470

第三版

免疫原:Fusion Protein of EIF2A

亚型:IgG1

克隆:1G7H4

描述:EIF2A,also named as CDA02,MSTP004 and MSTP089,belongs to the WD repeat EIF2A family. It functions in the early steps of protein synthesis of a small number of specific mRNAs. EIF2A,an alternative initiator tRNA-binding protein,has a key role in the translation of HCV mRNA during HCV infection,in turn promoting eIF2a phosphorylation by activating the eIF2a kinase PKR.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF