

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse Factor XIIIa Monoclonal Antibody

单克隆抗体

产品货号: MA02474

第三版

免疫原:Fusion protein of Factor XIIIa

亚型:IgG2b

克隆:2D2B2

**描述:**Factor XIIIa is a blood proenzyme that catalyzes intermolecular cross-linking of fibrinogen and is involved in clotting cascade. It has been identified in various tissues including placenta, uterus, and prostate. Anti- Factor XIIIa has been found to be useful in differentiating between dermatofibroma (90% positive) and desmoplastic malignant melanoma (negative). Factor XIIIa positivity is also seen in capillary hemangioblastoma, hemangioendothelioma, hemangiopericytoma, xanthogranuloma, xanthoma, hepatocellular carcinoma, glomus tumor, and meningioma.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF