



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat GAMT Monoclonal Antibody

单克隆抗体

产品货号: MA02478

第三版

免疫原:Fusion Protein of GAMT

亚型:IgG2a

克隆:1E7B11

描述:GAMT converts the latter to creatine with S-adenosylmethionine as the methyl donor. GAMT was induced in human and mouse cells in response to genotoxic and metabolic stress. GAMT was required for p53-dependent upregulation of creatine biosynthesis,fatty acid oxidation,and apoptosis in response to glucose deprivation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF