

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat HPX Monoclonal Antibody

单克隆抗体

产品货号: MA02485

第三版

免疫原:Fusion protein of Hemopexin

亚型:IgG1

克隆:3A9D6

描述:Hemopexin (HPX) is the plasma protein responsible for scavenging heme,thus preventing heme-mediated oxidative stress and heme-bound iron loss. In addition,hemopexin blocks heme activation of immune receptors and vascular inflammatory processes. It is mainly expressed in liver,the synthesis of which is induced after inflammation. Alterations of plasma hemopexin level have been linked to disorders like atherosclerosis and inflammatory diseases.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF