



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL21R Monoclonal Antibody

单克隆抗体

产品货号: MA02493

第三版

免疫原:Fusion Protein of IL21R

亚型:IgG2b

克隆:1D1C2

描述:IL21R,a receptor for interleukin-21,belongs to the type I cytokine receptor family. IL21R has been shown to form a heterodimeric receptor complex with the common gamma-chain,a receptor subunit also shared by the receptors for interleukin 2,4,7,9,and 15. This receptor transduces the growth promoting signal of IL21,and is important for the proliferation and differentiation of T cells,B cells,and natural killer (NK) cells. The ligand binding of this receptor leads to the activation of multiple downstream signaling molecules,including JAK1,JAK3,STAT1,and STAT3.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF