



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat LAP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02498

第三版

免疫原:Fusion Protein of LAP3

亚型:IgG1

克隆:1B11B8

**描述:**LAP3(Leucine aminopeptidase 3) is also named as LAPEP,PEPS(peptidase S) and belongs to the peptidase M17 family. It can catalyzes the removal of unsubstituted N-terminal amino-acids from various peptides and be presumably involved in the processing end regular turnover of intracellular proteins. PEPS is found in a wide variety of tissues and cultured cells but not in red cells and skin. The enzyme can utilize several di-,tri-,and tetrapeptides as substrates and is the slowest migrating of the peptidases after electrophoresis.. LAP3 has two isoform with MW of 52 and 56 kDa.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF