

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NFE2 Monoclonal Antibody

单克隆抗体

产品货号: MA02514

第三版

免疫原:Fusion Protein of NFE2

亚型:IgG1

克隆:1D7C8

描述:The transcription factor p45 nuclear factor-erythroid-derived 2(NFE2) is one component of the NF-E2 complex that essential for regulating erythroid and megakaryocytic maturation and differentiation. It also plays crucial roles in erythroid and megakaryocytic lineages. NFE2 recognizes the TCAT/C sequence of the AP-1-like core palindrome present in a number of erythroid and megakaryocytic gene promoters, and acts as an activator. And it may have a role in all aspects of hemoglobin production from globin and heme synthesis to procurement of iron.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF