



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat P4HB Monoclonal Antibody

单克隆抗体

产品货号: MA02528

第三版

免疫原:Fusion Protein of PDI

亚型:IgG2b

克隆:2E6A11

描述:PDIA1(Protein disulfide-isomerase) is also named as ERBA2L,PDI,P4HB,PO4DB. It is a multifunctional protein that catalyzes the formation,breakage and rearrangement of disulfide bonds. In some cell types,it seems to be secreted or associated with the plasma membrane,where it undergoes constant shedding and replacement from intracellular sources.It can exist as homodimer and monomers and homotetramers may also occur.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF