



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PGRMC1 Monoclonal Antibody

单克隆抗体

产品货号: MA02533

第三版

免疫原:Fusion Protein of PGRMC1

亚型:IgG1

克隆:1C10D11

描述:Progesterone receptor membrane component 1 (PGRMC1) is a member of a multi-protein progesterone-binding complex. However,PGRMC1 shares homology with cytochrome b5-related proteins rather than hormone receptors . It is a heme binding protein with biding sites for Src homology (SH2) and SH3 domain-containing proteins . PGRMC1 is overexpressed in a variety of cancers,and thus represents an important biomarker for cancer progression and a potential target for anticancer drugs . In nonmalignant tissues,PGRMC1 is highly expressed in the liver and kidney .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF