

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PLOD2 Monoclonal Antibody

单克隆抗体

产品货号: MA02541

第三版

免疫原:Fusion Protein of PLOD2

亚型:IgA

克隆:1H9E1

描述:PLOD2,also named as LH2,forms hydroxylysine residues in -Xaa-Lys-Gly- sequences in collagens. It is a potential novel prognostic factor for HCC patients following surgery. Among the PLOD genes,PLOD2 contributes to cancer prognosis and angiogenesis. Several authors have reported that PLOD2 expression might provide prognostic information about malignant tumours such as glioblastoma. PLOD2 expression is a useful biomarker for the effects of antiangiogenic treatment for malignancy.. It has 2 isoforms produced by alternative splicing and seven glycosylation sites.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF