



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PSMB10 Monoclonal Antibody

单克隆抗体

产品货号: MA02545

第三版

免疫原:Fusion Protein of PSMB10

亚型:IgG2a

克隆:1G8C1

描述:SMB10(Proteasome subunit beta type-10) is also named as LMP10,MECL1.It belongs to the peptidase T1B family and is a multicatalytic proteinase complex which is characterized by its ability to cleave peptides with Arg,Phe,Tyr,Leu,and Glu adjacent to the leaving group at neutral or slightly basic pH.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF