

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat RGS4 Monoclonal Antibody

单克隆抗体

产品货号: MA02547

第三版

免疫原:Fusion Protein of RGS4

亚型:IgG1

克隆:1E2E2

描述:Regulator of G protein signaling (RGS) proteins are a large family of highly diverse,multifunctional signaling proteins. RGS4 is one of seven members of a classic R4 RGS protein that accelerates the GTPase activity of the Gai/o and Gaq/11 family members. RGS4 is highly expressed in multiple brain regions including the cerebral cortex,amygdala,hippocampus,and striatum. Extensive studies have demonstrated that RGS4 plays an important role in regulating smooth muscle contraction,cardiomyocyte development,neural plasticity and psychiatric disorders.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF