

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SKAP2 Monoclonal Antibody

单克隆抗体

产品货号: MA02553

第三版

免疫原:Fusion Protein of SKAP2

亚型:IgG1

克隆:2G3G9

描述:SKAP2 (src kinase associated phosphoprotein 2),also named PRAP,contains PH domain and SH3 domain. It is localized in cytoplasm and was expressed in ubiquitously in various tissues including lung,skeletal muscle,and spleen,but not detectable in thymus and T cell lymphoma cells. SKAP2 may negatively regulate alpha-synuclein phosphorylation following cell stress. SKAP2 may be involved in B-cell and macrophage adhesion processes by coupling the B-cell receptor (BCR) to integrin activation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF