

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse SOX-2 Monoclonal Antibody

单克隆抗体

产品货号: MA02559

---

第三版

免疫原:Fusion Protein of SOX2

亚型:IgG1

克隆:4B9B10

描述:Sox2,also known as SRY (sex determining region Y)-box 2,is a transcription factor essential for maintaining self-renewal of undifferentiated ES cells and is one of the key transcription factors used to reprogram mouse and human fibroblasts to a pluripotent state. Sox2 expressed in undifferentiated pluripotent stem cells and germ cells during development.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF