



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SYCP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02565

第三版

免疫原:Fusion Protein of SYCP3

亚型:IgG1

克隆:1A11F1

描述:SYCP3 is an essential structural component of the synaptonemal complex. This complex is involved in synapsis, recombination and segregation of meiotic chromosomes. SYCP3 is required for centromere pairing during meiosis in male germ cells. SYCP3 is also required for normal meiosis during spermatogenesis and male fertility. Mutations in SYCP3 are associated with azoospermia in males and susceptibility to pregnancy loss in females.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF