



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TCN2 Monoclonal Antibody

单克隆抗体

产品货号: MA02568

第三版

免疫原:Fusion Protein of TCN2

亚型:IgG2b

克隆:2G6D2

描述:TCN2 (Transcobalamin-2) is a vitamin B12 (cobalamin; Cbl) binding and transport protein in plasma that facilitates the cellular uptake of the vitamin. TCN2 is secreted in a constitutive fashion in most cells. Plasma TCN2 levels have been reported to be elevated in a number of diverse human diseases. Defects in TCN2 are cause of transcobalamin-2 deficiency. (10714245)

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF