



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat TFF2 Monoclonal Antibody

单克隆抗体

产品货号: MA02569

第三版

免疫原:Fusion Protein of TFF2

亚型:IgG2a

克隆:2D7B1

描述:The trefoil factor family (TFF),comprises of three polypeptides,TFF1,TFF2 and TFF3 (7-12 kDa), secreted to mucosal surfaces by mucus producing cells,prominently in the gastrointestinal tract. TFF2,also known as spasmolytic polypeptide,is a low-molecular weight protein,expressed in mucous neck cells of the fundus and glands at the base of the antrum in normal human stomach. TFF2 could Inhibit gastrointestinal motility and gastric acid secretion. However,recent studies suggest that TFF2 could also play an important role in the immune system.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF