



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NSA2 Monoclonal Antibody

单克隆抗体

产品货号: MA02573

第三版

免疫原:Fusion Protein of TINP1

亚型:IgG1

克隆:1E3B11

描述:TINP1 (TGF beta inducible nuclear protein 1) was widely expressed in various cell lines. TINP1 overexpression significantly promoted tumor cell proliferation, which may be associated with the downregulation of p53 expression. Furthermore, TINP1 promoted a number of cell lines to the S phase of the cell cycle. TINP1 promotes cell proliferation and significantly reduces p53 and p21 expression.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF