



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse TRIM24 Monoclonal Antibody

单克隆抗体

产品货号: MA02577

第三版

免疫原:Fusion Protein of TRIM24

亚型:IgG1

克隆:1D6D9

描述:TRIM24,also named as RNF82,TIF1 and TIF1A,interacts selectively in vitro with the AF2-activating domain of the estrogen receptors. Association with DNA-bound estrogen receptors,it requires the presence of estradiol. It is a regulation protein of P53. Defects in TRIM24 are a cause of thyroid papillary carcinoma (TPC) wich is a common tumor of the thyroid that typically arises as an irregular,solid or cystic mass from otherwise normal thyroid tissue.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF