



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat VAMP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02580

第三版

免疫原:Fusion Protein of VAMP3/Cellubrevin

亚型:IgG2b

克隆:1D10D10

描述:VAMP3,also named as cellubrevin or synaptobrevin 3,is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. VAMP3 is a v-SNARE (soluble NSF-attachment protein receptor). Characterized by a common sequence called the SNARE motif,SNARE proteins are involved in membrane fusion and vesicular transport . VAMP3 resides in recycling endosomes and endosome-derived transport vesicles,involved in regulating membrane traffic. It is implicated in recycling of transferrin receptors,secretion of alpha-granules in platelets,and membrane trafficking during cell migration.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF