



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human VPS37A Monoclonal Antibody

单克隆抗体

产品货号: MA02583

第三版

免疫原:Fusion Protein of VPS37A

亚型:IgG1

克隆:1F9A4

描述:VPS37A,also named as HCRP1,belongs to the VPS37 family. It is a regulator of vesicular trafficking process. VPS37A is a member of the endosomal sorting complex required for transport (ESCRT) system,in upper motor neuron disease. ESCRT has been shown to play a central role in intracellular trafficking,in the maturation of multivesicular bodies and the sorting of ubiquitinated membrane proteins into internal luminal vesicles.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF