



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human VTN Monoclonal Antibody

单克隆抗体

产品货号: MA02586

第三版

免疫原:Fusion Protein of VTN

亚型:IgG1

克隆:1D2C8

描述:Vitronectin is a glycoprotein present in blood and in the extracellular matrix . Vitronectin is a multifunctional glycoprotein that mediates cell-to-substrate adhesion,inhibits the cytolytic action of the terminal complement cascade in vitro and binds to several serine protease inhibitors of the serpin family .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF