

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat COL4A1 Monoclonal Antibody

单克隆抗体

产品货号: MA02646

---

第三版

免疫原:Synthetic Peptide of Collagen IV

亚型:IgG

克隆:6F6

描述:This gene encodes the major type IV alpha collagen chain of basement membranes. Like the other members of the type IV collagen gene family, this gene is organized in a head-to-head conformation with another type IV collagen gene so that each gene pair shares a common promoter.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF