



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Phospho-STAT3 (Tyr705) Monoclonal Antibody

单克隆抗体

产品货号: MA02707

第三版

免疫原:Synthetic Peptide of Phospho-Stat3 (Y705)

亚型:IgG

克隆:1E9

描述:Signal transducer and transcription activator that mediates cellular responses to interleukins, KITLG/SCF, LEP and other growth factors. Once activated, recruits coactivators, such as NCOA1 or MED1, to the promoter region of the target gene (PubMed:17344214).

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF