
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

2, 6-Sialyltransferase (E41) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01004

第三版

描述:The histopathological picture of Alzheimer ' s disease is characterized by senile plaques and neurofibrillary tangles, and because the senile plaques form first, they are considered the initial lesion. Senile plaques are known to be formed by accumulation of β -amyloid peptide (A β).

抗原:Synthetic peptide for a part of Rat 2, 6-sialyltransferase (E41 Form).

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Reacts with Rat 2, 6-sialyltransferase (E41 Form). Not cross-react with Human or Mouse.