

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

14-3-3 Protein Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01008

第三版

描述:The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors.

抗原:Synthetic peptide for a part of Human 14-3-3 gamma.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:React to Human 14-3-3 ， ， ， ， ， 。