

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Adiponectin Receptor 1 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01010

第三版

描述:Adiponectin has been shown to be one of the adipokines possessing anti-diabetic, anti-atherosclerotic and anti-inflammatory actions. It was recently shown that the blood adiponectin levels were reduced in cases of obesity.

抗原:Synthetic peptide of a part of human Adiponectin Receptor 1.

配方:Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Not cross-react with Adiponectin Receptor 2. Cross-reacts with Mouse.