

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

14-3-3 Protein Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01016

---

第三版

**描述:**The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors.

**抗原:**Synthetic peptide for C-terminal of Human 14-3-3 .

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**Human 14-3-3 gamma specific, Non-crossreact to Human 14-3-3 , , , .