

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

14-3-3 Protein Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01020

第三版

描述:The 14-3-3 proteins are a family of conserved regulatory molecules expressed in all eukaryotic cells. A striking feature of the 14-3-3 proteins is their ability to bind a multitude of functionally diverse signaling proteins, including kinases, phosphatases, and transmembrane receptors.

抗原:Synthetic peptide for a part of N-terminal of Human 14-3-3 (MDKNELVQK).

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Human 14-3-3 zeta specific, Non-crossreact to Human 14-3-3 , , , .