
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AIMP2 (Aminoacyl tRNA synthetase complex-interacting multifunctional protein 2)

多克隆抗体

产品货号: PA01021

第三版

描述: AIMP2 is an important component for association of ARS multi-complex. AIMP2 is a potent tumor suppressor. AIMP2 inhibit c-myc activity in TGFbeta signaling, link c-IAP1 to TRAF2 by TNF-alpha stimulation and activate p53 under DNA-damaged condition.

抗原: Denatured full length native recombinant human AIMP2.

配方: Anti-AIMP2 rabbit serum containing 0.05% sodium azide for preservative.

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: This antibody reacts specifically with 38kDa AIMP2(JTV-1) protein.