
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

14-3-3 Protein (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01028

第三版

描述:The cell cycle checkpoint plays an important role in maintaining the integrity of cells. The G1 arrest is due to p53-mediated induction of the cyclin-dependent kinase inhibitor p21WAF1/CIP1/SD11, but the basis for the G2 arrest is unknown.

抗原:Synthetic peptide of c-terminus portion for human 14-3-3 Sigma protein.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Human 14-3-3 Sigma specific. Non-cross react with 14-3-3 Beta, Gamma, Epsilon, Zeta, Eta, and Tau. (Confirmed by western blotting).