

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Aminoacyl tRNA synthetase complex interacting multifunctional protein 1 (AIMP1)

多克隆抗体

产品货号: PA01040

---

第三版

**描述:**AIMP1, previously known as p43, is an auxiliary component of ARS multi complex and a precursor of cytokine endothelial-monocyte activating polypeptide II (EMAPII). AIMP1 is known as a cytokine working in the control of angiogenesis, inflammation and wound healing. AIMP1 also has glucagon-like hormonal activity.

**抗原:**Native full length native recombinant human AIMP1.

**配方:**Anti-AIMP1 rabbit serum containing 0.05% sodium azide for preservative.

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**This antibody reacts specifically with 38kDa p43(AIMP1) protein.