

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Amyloid Beta ( N ) A Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01047

第三版

**描述:**Alzheimer ' s Disease (AD) is characterized by the presence of extracellular plaques and intracellular neurofibrillary tangles (NFTs) in the brain. The major protein component of these plaques is beta amyloid peptide (A<sub>β</sub>), a 40 to 43 amino acid peptide cleaved from amyloid precursor protein by  $\beta$ -secretase and a putative  $\gamma$  secretase.

**抗原:**Synthetic peptide for N-terminal of Human Amyloid Beta A<sub>β</sub> 1-42.

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**React with all of Human Amyloid  $\beta$  [1-40], [1-42] and [1-43].