
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Angiotensinogen (405) Anti-Mouse/Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01050

第三版

描述:Angiotensinogen is the precursor of angiotensin and is cleaved into angiotensin I and II in the renin-angiotensin system, and it has long been reported to play an important role in controlling blood pressure.

抗原:Synthetic peptide of the C terminal part of Mouse Angiotensinogen (IGDTNPRVGEVLNSIL).

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Reacts with mouse and rat Angiotensinogen.