

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

ApoE (A299) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01052

---

第三版

**描述:** Apolipoprotein-E (Apo-E) is a structural component of very low-density lipoprotein (VLDL) synthesized by the liver and intestinally synthesized chylomicrons.

**抗原:** Synthetic peptide for c-terminal of Human ApoE.

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Can be react to ApoE2, ApoE3 and ApoE4.