
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

APP (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01058

第三版

描述: Amyloid precursor protein (APP) is precursor protein of Amyloid which is major constituent of senile plaque in Alzheimer ' s disease. It is known that there are three major isoforms, APP695, APP751 and APP770, and are generated from alternative splicing of common precursor mRNA.

抗原: Synthetic peptide in portion of C-terminus of Human APP.

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃.

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Cross-reacts with Mouse and Rat. All isoforms (APP695, APP751, APP770) are detectable. Non-crossreacts to sAPP and sAPP .