

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

sAPP -Wild Type Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01064

第三版

**描述:** Amyloid precursor protein (APP) is precursor protein of Amyloid which is major constituent of senile plaque in Alzheimer ' s disease. It is known that there are three major isoforms, APP695, APP751 and APP770, and are generated from alternative splicing of common precursor mRNA.

**抗原:** Synthetic peptide in portion of C-terminus of Human sAPP -Wild Type (ISEVKM).

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>.

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** This antibody can detect sAPP-Wild type which is cleaved by -secretase. Non-crossreacts with sAPP and full-length APP.