

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-BACE1 (42) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01065

第三版

**描述:**The current understanding of the role of amyloid in Alzheimer's disease (AD) has been established by a remarkable congruence of multiple disciplines: neuropathology, biochemistry, molecular biology, and epidemiological genetics.

**抗原:**Synthetic peptide of a part of human BACE1 (RLPRETDEEPEEPGRR)

**配方:**Lyophilized product from 1 % BSA in PBS containing 0.05 % Na<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Reacts with human, mouse and rat BACE1.