



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Big Endothelin-1 (22-38) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01069

第三版

描述:Endothelins are 21 amino acid peptides that are produced ubiquitously by vascular endothelial cells, smooth muscle cells, and other cells in different organs. Endothelins are secreted as big-endothelins that are converted to active proteins by the endothelin-converting enzyme.

抗原:Synthetic peptides of rat big endothelin-1 (22-38)

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: