
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

BNP (Brain Natriuretic Peptide) Anti-Human Rabbit IgG

多克隆抗体

产品货号: PA01070

第三版

描述: This antibody can be stained in bouin fixed paraffin embedded tissues by several Immunohistochemical techniques. The optimal dilution is 4 – 20 ug/mL , however, the dilution rate should be optimized by each laboratories.

抗原: Synthetic Human BNP

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Confirmed in Human Brain