

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

C-ERC/Mesothelin (306) Anti-Rat Rabbit Igg Affinity Purify

多克隆抗体

产品货号: PA01071

---

第三版

**描述:** Erc has been identified as a gene showing stronger expression in cancer-affected renal cell than in normal renal tissue in Eker rats (a rat model of renal cancer). The human homologue of the protein encoded by this gene is called MPF (megakaryocyte potentiating factor) or mesothelin.

**抗原:** Synthetic peptide of the part of Rat C-ERC/Mesothelin (GKEPNVV DENLIFYQNWELEA)

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**