

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-c-Kit (K963) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01075

第三版

描述:The proto-oncogene c-kit encodes a transmembrane tyrosine kinase receptor, and its ligand for c-kit receptor has been identified as the stem cell factor (SCF). Recent experimental studies have shown that c-kit plays a key role in the development of a component of the pacemaker system that is required for generation of autonomic gut motility.

抗原:Synthetic peptide for C-terminal of Human c-Kit

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Human c-Kit (M.W.145kDa) specific