

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-c-Met (Y1235 Phosphorylated) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01076

第三版

描述:Overexpression of the hepatocyte growth factor receptor (c-Met/HGF receptor), a transmembrane tyrosine kinase encoded by the MET proto-oncogene, is involved in transformation and invasive behavior of human carcinomas and sarcomas.

抗原:Synthetic peptide of the phosphorylated Tyr1235 part of human c-Met. (YDKEY(pY)SVHNK)

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: