

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

c-Met Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01078

第三版

**描述:**Overexpression of the hepatocyte growth factor receptor (c-Met/HGF receptor), a transmembrane tyrosine kinase encoded by the MET proto-oncogene, is involved in transformation and invasive behavior of human carcinomas and sarcomas.

**抗原:**Synthetic peptide for C-terminal of Human c-Met

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**React to c-Met beta chain.